

CD300A

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD300A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10006

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11314

Isotype: IgG

Swiss Prot: Q9UGN4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the CD300 glycoprotein family of cell surface proteins found on leukocytes involved in immune response signaling pathways. This gene is located on chromosome 17 in a cluster with all but one of the other family members. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]

To place an order, please [Click HERE](#).