

NEDD8

Reactivity: Human

Tested applications: WB IHC ICC IF FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50- 1:200 ICC 1:50- 1:200 IF 1:20- 1:100
FC 1:20- 1:50

Calculated MW: 9kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NEDD8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NEDD-8

Catalog #: A10008

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 4738

Isotype: IgG

Swiss Prot: Q15843

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).