

NPC1L1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 149kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NPC1L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10049

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29881

Isotype: IgG

Swiss Prot: Q9UHC9

Purity: Affinity purification

For research use only.