

CR1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 224kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10059

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1378

Isotype: IgG

Swiss Prot: P17927

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).