

IL17RB

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL17RB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10147

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55540

Isotype: IgG

Swiss Prot: Q9NRM6

Purity: Affinity purification

For research use only.