

KLHL9

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:69kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KLHL9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A10149

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55958

Isotype:IgG

Swiss Prot:Q9P2J3

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).