

MMP1

Reactivity:Human

Tested applications:WB IHC ICC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human MMP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10168

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4312

Isotype:IgG

Swiss Prot:P03956

Purity:Affinity purification

For research use only.