

PDX1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human PDX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10173

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3651

Isotype: IgG

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).