

HLA-G

Reactivity: Human Mouse Rat

Tested applications: WB

Calculated MW: 338

Observed MW: Refer to figures

Immunogen:

Recombinant Protein

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10182

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 3135

Isotype: IgG

Swiss Prot: P17693

Purity: Affinity purification

For research use only.