

FKRP

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:3000

Calculated MW:55kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FKRP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A10194

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79147

Isotype:IgG

Swiss Prot:Q9H9S5

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein which is targeted to the medial Golgi apparatus and is necessary for posttranslational modification of dystroglycan. Mutations in this gene have been associated with congenital muscular dystrophy, mental retardation, and cerebellar cysts. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Oct 2008]

To place an order, please [Click HERE](#).