

SOX17

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50- 1:200 IF 1:20- 1:100

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SOX17

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10198

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64321

Isotype: IgG

Swiss Prot: Q9H6I2

Purity: Affinity purification

For research use only.