

CDH10

Reactivity:Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:88kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CDH10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A10204

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1008

Isotype:IgG

Swiss Prot:Q9Y6N8

Purity:Affinity purification

For research use only.

Background:

This gene encodes a type II classical cadherin of the cadherin superfamily. Alternative splicing of this gene results in multiple transcript variants. At least one of these variants encodes a preproprotein that is proteolytically processed to generate the mature cadherin protein. These integral membrane proteins mediate calcium-dependent cell-cell adhesion and are composed of a large N-terminal extracellular domain, a single membrane-spanning domain, and a small, highly conserved C-terminal cytoplasmic domain. The extracellular domain consists of 5 subdomains, each containing a cadherin motif, and appears to determine the specificity of the protein's homophilic cell adhesion activity. Type II (atypical) cadherins are defined based on their lack of a histidine-alanine-valine (HAV) cell adhesion recognition sequence specific to type I cadherins. This particular cadherin is predominantly expressed in brain and is putatively involved in synaptic adhesions, axon outgrowth and guidance. Mutations in this gene may be associated with lung squamous cell carcinoma and colorectal cancer in human patients.

To place an order, please [Click HERE](#).