

GOLGA4

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 261kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GOLGA4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GCP2; GOLG; p230; CRPF46; MU-RMS-40.18;

Catalog #: A10216

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2803

Isotype: IgG

Swiss Prot: Q13439

Purity: Affinity purification

For research use only.

Background:

The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. This gene encodes one of the golgins, a family of proteins localized to the Golgi. This protein has been postulated to play a role in Rab6-regulated membrane-tethering events in the Golgi apparatus. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.

To place an order, please [Click HERE](#).