

LY6E

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 14kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LY6E

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RIGE; SCA2; RIG-E; SCA-2; TSA-1;

To place an order, please [Click HERE](#).

Catalog #: A10225

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4061

Isotype: IgG

Swiss Prot: Q16553

Purity: Affinity purification

For research use only.