

TMPRSS15

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 113kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TMPRSS15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ENTK; PRSS7;

Catalog #: A10239

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5651

Isotype: IgG

Swiss Prot: P98073

Purity: Affinity purification

For research use only.

Background:

This gene encodes an enzyme that converts the pancreatic proenzyme trypsinogen to trypsin, which activates other proenzymes including chymotrypsinogen and procarboxypeptidases. The precursor protein is cleaved into two chains that form a heterodimer linked by a disulfide bond. This protein is a member of the trypsin family of peptidases. Mutations in this gene cause enterokinase deficiency, a malabsorption disorder characterized by diarrhea and failure to thrive.

To place an order, please [Click HERE](#).