

RORC

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RORC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOR; RORG; RZRG; IMD42; NR1F3; RZR-GAMMA;

Catalog #: A10240

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6097

Isotype: IgG

Swiss Prot: P51449

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a DNA-binding transcription factor and is a member of the NR1 subfamily of nuclear hormone receptors. The specific functions of this protein are not known; however, studies of a similar gene in mice have shown that this gene may be essential for lymphoid organogenesis and may play an important regulatory role in thymopoiesis. In addition, studies in mice suggest that the protein encoded by this gene may inhibit the expression of Fas ligand and IL2. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).