

SMARCA1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 123kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SMARCA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SWI; ISWI; SWI2; SNF2L; SNF2L1; SNF2LB; SNF2LT; hSNF2L; NURF140;

Catalog #: A10248

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6594

Isotype: IgG

Swiss Prot: P28370

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SWI/SNF family of proteins. The encoded protein is an ATPase which is expressed in diverse tissues and contributes to the chromatin remodeling complex that is involved in transcription. The protein may also play a role in DNA damage, growth inhibition and apoptosis of cancer cells. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).