

TSPAN4

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TSPAN4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NAG2; NAG-2; TM4SF7; TSPAN-4; TETRASPAN;

Catalog #: A10253

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7106

Isotype: IgG

Swiss Prot: O14817

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein and is similar in sequence to its family member CD53 antigen. It is known to complex with integrins and other transmembrane 4 superfamily proteins. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).