

XPNPEP2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human XPNPEP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APP2; AEACEI;

Catalog #: A10255

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7512

Isotype: IgG

Swiss Prot: O43895

Purity: Affinity purification

For research use only.

Background:

Aminopeptidase P is a hydrolase specific for N-terminal imido bonds, which are common to several collagen degradation products, neuropeptides, vasoactive peptides, and cytokines. Structurally, the enzyme is a member of the 'pita bread fold' family and occurs in mammalian tissues in both soluble and GPI-anchored membrane-bound forms. A membrane-bound and soluble form of this enzyme have been identified as products of two separate genes.

To place an order, please [Click HERE](#).