

EIF3B

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 92kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EIF3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRT1; EIF3S9; EIF3-ETA; EIF3-P110; EIF3-P116;

To place an order, please [Click HERE](#).

Catalog #: A10259

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8662

Isotype: IgG

Swiss Prot: P55884

Purity: Affinity purification

For research use only.