

USP13

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 97kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human USP13

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ISOT3; IsoT-3;

To place an order, please [Click HERE](#).

Catalog #: A10264

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8975

Isotype: IgG

Swiss Prot: Q92995

Purity: Affinity purification

For research use only.