

SLC23A1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 65kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC23A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SVCT1; YSPL3; SLC23A2;

Catalog #: A10276

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9963

Isotype: IgG

Swiss Prot: Q9UHI7

Purity: Affinity purification

For research use only.

Background:

The absorption of vitamin C into the body and its distribution to organs requires two sodium-dependent vitamin C transporters. This gene encodes one of the two transporters. The encoded protein is active in bulk vitamin C transport involving epithelial surfaces. Previously, this gene had an official symbol of SLC23A2. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).