

FSTL3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FSTL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FLRG; FSRP;

Catalog #: A10279

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10272

Isotype: IgG

Swiss Prot: O95633

Purity: Affinity purification

For research use only.

Background:

Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis.

To place an order, please [Click HERE](#).