

CD226

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD226

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PTA1; DNAM1; DNAM-1; TLISA1;

Catalog #: A10282

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10666

Isotype: IgG

Swiss Prot: Q15762

Purity: Affinity purification

For research use only.

Background:

This gene encodes a glycoprotein expressed on the surface of NK cells, platelets, monocytes and a subset of T cells. It is a member of the Ig-superfamily containing 2 Ig-like domains of the V-set.

The protein mediates cellular adhesion of platelets and megakaryocytic cells to vascular endothelial cells. The protein also plays a role in megakaryocytic cell maturation. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).