

FNBP1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FNBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FBP17;

Catalog #: A10289

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23048

Isotype: IgG

Swiss Prot: Q96RU3

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the formin-binding-protein family. The protein contains an N-terminal Fer/Cdc42-interacting protein 4 (CIP4) homology (FCH) domain followed by a coiled-coil domain, a proline-rich motif, a second coiled-coil domain, a Rho family protein-binding domain (RBD), and a C-terminal SH3 domain. This protein binds sorting nexin 2 (SNX2), tankyrase (TNKS), and dynamin; an interaction between this protein and formin has not been demonstrated yet in human.

To place an order, please [Click HERE](#).