

## ICMT

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 32kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human ICMT

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PCMT; PPMT; PCCMT; HSTE14; MST098; MSTP098;

**Catalog #:** A10293

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 23463

**Isotype:** IgG

**Swiss Prot:** O60725

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes the third of three enzymes that posttranslationally modify isoprenylated C-terminal cysteine residues in certain proteins and target those proteins to the cell membrane.

This enzyme localizes to the endoplasmic reticulum. Alternative splicing may result in other transcript variants, but the biological validity of those transcripts has not been determined.

*To place an order, please [Click HERE](#).*