

TRHDE

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:117kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TRHDE

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAP-II; PGPEP2; TRH-DE;

Catalog #:A10300

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:29953

Isotype:IgG

Swiss Prot:Q9UKU6

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the peptidase M1 family. The encoded protein is an extracellular peptidase that specifically cleaves and inactivates the neuropeptide thyrotropin-releasing hormone.

To place an order, please [Click HERE](#).