

IL20RA

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 62kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL20RA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRF2-8; IL-20R1; IL-20RA; IL-20R-alpha;

Catalog #: A10308

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 53832

Isotype: IgG

Swiss Prot: Q9UHF4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the type II cytokine receptor family. The encoded protein is a subunit of the receptor for interleukin 20, a cytokine that may be involved in epidermal function. The interleukin 20 receptor is a heterodimeric complex consisting of the encoded protein and interleukin 20 receptor beta. This gene and interleukin 20 receptor beta are highly expressed in skin, and are upregulated in psoriasis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).