

GNG2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 8kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GNG2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10309

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54331

Isotype: IgG

Swiss Prot: P59768

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of the gamma subunits of a guanine nucleotide-binding protein. Such proteins are involved in signaling mechanisms across membranes. Various subunits forms heterodimers which then interact with the different signal molecules.

To place an order, please [Click HERE](#).