

NLGN3

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:94kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NLGN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNL3;

Catalog #:A10310

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54413

Isotype:IgG

Swiss Prot:Q9NZ94

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. Mutations in this gene may be associated with autism and Asperger syndrome. Multiple transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).