

GPR126

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:137kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GPR126

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APG1; DREG; VIGR; LCCS9; GPR126; PS1TP2;

Catalog #:A10318

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57211

Isotype:IgG

Swiss Prot:Q86SQ4

Purity:Affinity purification

For research use only.

Background:

This gene, which is upregulated in human umbilical vein endothelial cells, encodes a G protein-coupled receptor. Variations in this gene can affect a person's stature. Multiple transcript variants encoding different proteins have been found for this gene.

To place an order, please [Click HERE](#).