

LRRC4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 73kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LRRC4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NAG14; NGL-2;

Catalog #: A10321

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64101

Isotype: IgG

Swiss Prot: Q9HBW1

Purity: Affinity purification

For research use only.

Background:

This gene is significantly downregulated in primary brain tumors. The exact function of the protein encoded by this gene is unknown.

To place an order, please [Click HERE](#).