

SPG11

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 279kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SPG11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALS5; CMT2X; KIAA1840;

Catalog #: A10326

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80208

Isotype: IgG

Swiss Prot: Q96J17

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a potential transmembrane protein that is phosphorylated upon DNA damage. Defects in this gene are a cause of spastic paraplegia type 11 (SPG11). Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).