

FCRL4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FCRL4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FCRH4; IGFP2; IRTA1; CD307d;

Catalog #: A10329

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83417

Isotype: IgG

Swiss Prot: Q96PJ5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the immunoglobulin receptor superfamily and is one of several Fc receptor-like glycoproteins clustered on the long arm of chromosome 1. The encoded protein has four extracellular C2-type immunoglobulin domains, a transmembrane domain and a cytoplasmic domain that contains three immune-receptor tyrosine-based inhibitory motifs. This protein may play a role in the function of memory B-cells in the epithelia. Aberrations in the chromosomal region encoding this gene are associated with non-Hodgkin lymphoma and multiple myeloma.

To place an order, please [Click HERE](#).