

RAB34

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 29kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RAB34

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAH; RAB39

Catalog #: A10332

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83871

Isotype: IgG

Swiss Prot: P0DI83

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein belonging to the RAB family of proteins, which are small GTPases involved in protein transport. This family member is a Golgi-bound member of the secretory pathway that is involved in the repositioning of lysosomes and the activation of macropinocytosis. Alternative splicing of this gene results in multiple transcript variants. This gene overlaps and shares exon structure with the nine-amino acid residue-repeats (NARR) gene, which encodes a functionally distinct nucleolar protein from a different reading frame.

To place an order, please [Click HERE](#).