

CABLES1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:68kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CABLES1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CABL1; IK3-1; CABLES; HsT2563;

Catalog #:A10335

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:91768

Isotype:IgG

Swiss Prot:Q8TDN4

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein involved in regulation of the cell cycle through interactions with several cyclin-dependent kinases. One study (PMID: 16177568) reported aberrant splicing of transcripts from this gene which results in removal of the cyclin binding domain only in human cancer cells, and reduction in gene expression was shown in colorectal cancers (PMID: 17982127). Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).