

SNX18

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SNX18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SNAG1; SH3PX2; SH3PXD3B;

Catalog #: A10337

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 112574

Isotype: IgG

Swiss Prot: Q96RF0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a SH3 domain. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).