

DNAJB2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DNAJB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HSJ1; CMT2T; DSMA5; HSJ-1; HSPF3;

Catalog #:A10345

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3300

Isotype:IgG

Swiss Prot:P25686

Purity:Affinity purification

For research use only.

Background:

This gene is almost exclusively expressed in the brain, mainly in the neuronal layers. It encodes a protein that shows sequence similarity to bacterial DnaJ protein and the yeast homologs. In bacteria, this protein is implicated in protein folding and protein complex dissociation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).