

KIAA0101

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:12kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human KIAA0101

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

L5; PAF; OEATC; PAF15; OEATC1; p15PAF; NS5ATP9; OEATC-1; p15/PAF; p15(PAF)

To place an order, please [Click HERE](#).

Catalog #:A10357

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9768

Isotype:IgG

Swiss Prot:Q15004

Purity:Affinity purification

For research use only.