

SERPINA5

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 46kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic Peptide of human SERPINA5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SERPINA5;PAI3;PCI;PLANH3;PROCI ;

Catalog #: A1036

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5104

Isotype: IgG

Swiss Prot: P05154

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. This gene is one in a cluster of serpin genes located on the q arm of chromosome 14. This family member is a glycoprotein that can inhibit several serine proteases, including protein C and various plasminogen activators and kallikreins, and it thus plays diverse roles in hemostasis and thrombosis in multiple organs. [provided by RefSeq, Aug 2012]

To place an order, please [Click HERE](#).