

NOP14

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 97kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NOP14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UTP2; NOL14; C4orf9; RES425; RES4-25

Catalog #: A10361

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8602

Isotype: IgG

Swiss Prot: P78316

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that plays a role in pre-18s rRNA processing and small ribosomal subunit assembly. The encoded protein may be involved in the regulation of pancreatic cancer cell proliferation and migration. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).