

ADAM10

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 98kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ADAM10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAK; kuz; AD10; AD18; MADM; CD156c; HsT18717

Catalog #: A10438

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 102

Isotype: IgG

Swiss Prot: O14672

Purity: Affinity purification

For research use only.

Background:

Members of the ADAM family are cell surface proteins with a unique structure possessing both potential adhesion and protease domains. This gene encodes an ADAM family member that cleaves many proteins including TNF-alpha and E-cadherin.

To place an order, please [Click HERE](#).