

PRLHR

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:41kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PRLHR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GR3; GPR10; PrRPR;

Catalog #:A10479

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2834

Isotype:IgG

Swiss Prot:P49683

Purity:Affinity purification

For research use only.

Background:

PRLHR is a 7-transmembrane domain receptor for prolactin-releasing hormone (PRLH; MIM 602663) that is highly expressed in anterior pituitary (Ozawa et al., 2002

To place an order, please [Click HERE](#).