

INTS9

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 74kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human INTS9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

INT9; RC74; CPSF2L;

Catalog #: A10480

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55756

Isotype: IgG

Swiss Prot: Q9NV88

Purity: Affinity purification

For research use only.

Background:

This gene encodes a subunit of the Integrator complex. This protein complex binds the C-terminal domain of RNA polymerase II and likely plays a role in small nuclear RNA processing. The encoded protein has similarities to the subunits of the cleavage and polyadenylation specificity factor complex. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).