

CPB1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CPB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CPB; PASP; PCPB;

Catalog #: A10486

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1360

Isotype: IgG

Swiss Prot: P15086

Purity: Affinity purification

For research use only.

Background:

Three different procarboxypeptidases A and two different procarboxypeptidases B have been isolated. The B1 and B2 forms differ from each other mainly in isoelectric point. Carboxypeptidase B1 is a highly tissue-specific protein and is a useful serum marker for acute pancreatitis and dysfunction of pancreatic transplants. It is not elevated in pancreatic carcinoma.

To place an order, please [Click HERE](#).