

DNAJC5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DNAJC5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSP; NCL; CLN4; CLN4B; DNAJC5A; mir-941-2; mir-941-3; mir-941-4; mir-941-5;

Catalog #: A10489

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80331

Isotype: IgG

Swiss Prot: Q9H3Z4

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the J protein family. J proteins function in many cellular processes by regulating the ATPase activity of 70 kDa heat shock proteins. The encoded protein plays a role in membrane trafficking and protein folding, and has been shown to have anti-neurodegenerative properties. The encoded protein is known to play a role in cystic fibrosis and Huntington's disease. A pseudogene of this gene is located on the short arm of chromosome 8.

To place an order, please [Click HERE](#).