

MTAP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MTAP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MTAP;MSAP;c86fus ;

Catalog #: A1049

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4507

Isotype: IgG

Swiss Prot: Q13126

Purity: Affinity purification

For research use only.

Background:

MTAP is an enzyme that is essential for the salvage pathway for both adenine and methionine synthesis. MTAP catalyzes the cleavage of 5-methylthioadenosine into adenine and 5-methylthio-D-ribose-1-phosphate. Adenine is then used to generate AMP whereas 5-methylthio-D-ribose-1-phosphate is converted into methionine (1,2). MTAP is expressed in all normal cells and tissues, although frequently lost in different human tumors including pancreatic adenocarcinoma, neuroendocrine tumors, non-small cell lung carcinoma and breast carcinoma. MTAP is usually codeleted with p16 (cdkN2a/ARF) (3-5). MTAP overexpression in breast cancer cells inhibits their ability to form colonies in soft agar, thereby implicating its function as a tumor suppressor (6).

To place an order, please [Click HERE](#).