

ADAM8

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 89kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADAM8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MS2; CD156; CD156a;

Catalog #: A10497

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 101

Isotype: IgG

Swiss Prot: P78325

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene may be involved in cell adhesion during neurodegeneration, and it is thought to be a target for allergic respiratory diseases, including asthma. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).