

CLEC4G

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:33kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CLEC4G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LP2698; UNQ431; DTTR431; LSECTin;

Catalog #:A10500

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:339390

Isotype:IgG

Swiss Prot:Q6UXB4

Purity:Affinity purification

For research use only.

Background:

This gene encodes a glycan-binding receptor and member of the C-type lectin family which plays a role in the T-cell immune response. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).