

CSF2RB

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 97kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CSF2RB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD131; IL3RB; IL5RB; SMDP5; CDw131

Catalog #: A10501

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1439

Isotype: IgG

Swiss Prot: P32927

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is the common beta chain of the high affinity receptor for IL-3, IL-5 and CSF. Defects in this gene have been reported to be associated with protein alveolar proteinosis (PAP).

To place an order, please [Click HERE](#).